

# Oxy HOME<sup>®</sup>

## Stationary Concentrator

Keep Going.<sup>sm</sup>

NEW

### BENEFITS

- Quiet Operation and Small Footprint
- Purity Check by Percentage
- Alarm Indicators: No Flow, Low/High Pressure, Purity, Power
- Twin Blower Cooling Design
- Contemporary Design
- 3-Year Replacement Warranty

### PRODUCT SPECS **Part Number 1400-5100**

|                                |  |
|--------------------------------|--|
| <b>Flow Specifications</b>     | 0 - 5 LPM  |
| <b>Oxygen Concentration</b>    | 93% (+3%, -3%)                                       |
| <b>Oxygen Output Pressure</b>  | 8 (+2, -2) psi                                       |
| <b>Rated Voltage/Frequency</b> | 120 VAC ±10% 60Hz                                    |
| <b>Power Consumption</b>       | 350 Watts  |
| <b>Sound</b>                   | ≤ 45 dB(A)   |
| <b>Dimensions</b>              | 13.4" W x 11.8" D x 25.6" H<br>34 cm x 30 cm x 65 cm |
| <b>Weight</b>                  | 41.8 lbs (19 kg)                                     |



# 1400-5100

CONTACT US TODAY TO LEARN MORE

888-327-7301 • [orders@oxygo.life](mailto:orders@oxygo.life) • [oxygo.life](http://oxygo.life)